

Wednesday 6 May 2015

**ANTEO TECHNOLOGIES INTRODUCES NEW PRODUCT TO
MARKET
AMG™ Coupling Kit, 200 nm Magnetic Particles**

Today Anteo Technologies Pty Ltd (a subsidiary of Anteo Diagnostics Limited, ASX: ADO) has launched a new product - AMG™ Coupling Kit, 200 nm Magnetic Particles - to market.

The Anteo Mix&Go™ (AMG) kit has been specifically developed by Anteo's product development team to expand the AMG line of Mix&Go enabled products. The range offers novel advantages over traditional methods. This 200 nm Magnetic Particle kit gives researchers the unique ability to bind a range of biomolecules onto small magnetic particles. The smaller particles provide a number of benefits when used in immunoassay applications, including lateral flow assays, and when used for biomolecular separations including cell separation applications. This product will deliver scientists greater control and accuracy in less time.

Dr Geoff Cumming, CEO, Anteo Technologies said: "Our newest Coupling Kit effectively increases the surface area per unit mass available for protein binding. This can lead to a number of improvements including increased assay sensitivity. Normally, the handling complexity increases disproportionately with reducing size. Using Mix&Go, we have overcome the tendency for smaller sized particles to aggregate and fall out of solution. This is a unique feature of Mix&Go."

Anteo's patented Mix&Go technology platform, which is used across the company's suite of AMG products, is renowned for providing stability and excellent reproducibility of results. Anteo's platform technology allows activated particles to remain stable for at least one year.

How the AMG™ Coupling Kit, 200 nm Magnetic Particles works

The particles in the AMG™ Coupling Kit are provided fully activated with Mix&Go and are ready to use, with general buffer solutions known to work for a majority of applications.

The kit allows users to immediately couple with their protein of choice.

The kit is suitable for coupling a broad range of biomolecules for diverse applications, such as magnetic, fluorescent and chemiluminescent detection in immunoassays and purification applications using biomolecules including proteins, peptides, antigens and antibodies.

For personal use only

Anteo's suite of AMG products

The AMG Coupling Kit, 200 nm Magnetic Particles is the second coupling kit in the AMG line and the fifth product released to market in the range. Anteo has now delivered a number of products exemplifying (1) protein bound, (2) activation ready and (3) coupling ready products. Anteo is now focused on using the AMG product range to deliver on its commercialisation strategy.

Anteo Technologies

www.anteotech.com

Product information and online orders: www.anteotech.com/product/lnpckmp

Anteo Technologies Pty Ltd

Anteo uses its patented technology to develop, manufacture and commercialise proprietary surface coatings for use in healthcare, life sciences and beyond. Its patented technology is applied in the Mix&Go product range, which delivers solutions to the challenges of establishing highly functional biological interfaces between fragile biomolecules, and synthetic, and often incompatible, materials.

For personal use only